

90D



PRIMARY APPLICATION

**Pan Retinal Exam and
Small Pupil Examination**

- + Original 90D lens that started the slit lamp fundus examination revolution and the industry gold standard
- + Small diameter ring is ideal for dynamic fundoscopy and easy manipulation within the orbit
- + Optical profile makes it easy to use - an ideal training lens for new students and residents
- + Outstanding general diagnostic lens for pan retinal examination
- + Can be used on small pupils and patients who do not accommodate dilation

74°/89°
FIELD OF
VIEW

0.76x
IMAGE
MAG

1.32x
LASER
SPOT MAG

